BIO-C[®] REPAIR

CHOOSE YOUR LANGUAGE

ENGLISH.....	3	SLOVENŠČINA	65
ESPAÑOL.....	8	SLOVENSKY	70
PORTUGUÊS	14	POLSKI	76
FRANÇAIS	19	LIETUVIŠKAI	81
ITALIANO	25	ČESKY	87
ROMÂNĂ	31	MAGYAR	92
NEDERLANDS.....	36	TÜRKÇE	98
DEUTSCH	42	БЪЛГАРСКИ	103
SVENSKA	48	ΕΛΛΗΝΙΚΑ	109
SUOMI	54	日本	116
HRVATSKI	59		

ENGLISH

DEFINITION

Bioceramic reparative cement ready for use.

INDICATIONS

- A. Treatment of perforation of root or furcation by canal;
- B. Treatment of perforation of root or furcation surgically;
- C. Treatment of internal reabsorption by canal;
- D. Treatment of external or communicating internal reabsorption surgically;
- E. Retrofilling in paraendodontic surgery;
- F. Direct and indirect pulp capping;
- G. Apexification;
- H. Apicogenesis and Pulpotomy.
- I. Pulp regeneration

COMPOSITION

Calcium silicates, calcium aluminate, calcium oxide, zirconium oxide, iron oxide, silicon dioxide and dispersing agent.

TECHNICAL INFORMATION

- Bioactive product ready for use;
- Work Time: product ready for use, mixture not being necessary;
- Setting Time: \leq 120 minutes, depending upon the humidity of the place;
 - o Important: dry the cavity to be applied with paper cones, but without drying excessively; the humidity coming from the tissue will start the setting reaction of the material;
- High alkalinity - pH: \sim 12;
- High radiopacity: \geq 7 mm of the aluminum scale (standard ISO 6876);
- Threadable syringe to facilitate the procedure of removing the product from the syringe.

Important: Before use check if the product is hardened at the tip of the syringe and, if so, discard the first portion.

USE TECHNIQUE

A. TREATMENT OF PERFORATION OF ROOT OR FURCATION BY CANAL

1. Anesthetize, install the absolute isolation, execute the biomechanical preparation of the duct and cleaning

of the existing perforation;

2. Fill the canal in a conventional way leaving the perforation exposed;
3. Execute hemostatic control and maintain the cavity slightly moist;
4. Insert BIO-C® REPAIR and, if there is overflow, remove the excess material;
5. X-ray to check the correct seating of the material in the place applied;
6. Fill in the remainder of the duct with guttapercha cones and a filling cement such as MTA-FILLAPEX or BIO-C® SEALER;
7. Execute the coronary sealing and restoration.

Important:

- In perforations of furcation apply a layer of glass ionomer as a mechanical barrier followed by restoration with final material;
- Do not use composite resin directly on BIO-C REPAIR before its final setting.
- X-ray and execute preservation for at least two years.

B. TREATMENT OF PERFORATION OF ROOT OR FURCATION SURGICALLY

1. Anesthetize, raise the gingival flap and execute osteotomy to locate the perforation;
2. Prepare the perforation with an appropriate instrument to facilitate the insertion of the material;
3. Execute hemostatic control and maintain the cavity slightly moist;
4. Insert BIO-C® REPAIR in the cavity, adapt with condensers and remove excess if overflowing;
5. Reposition the gingival flap and sew it;
6. X-ray and execute preservation for at least two years;

C. TREATMENT OF INTERNAL REABSORPTION BY CANAL

1. Anesthetize and install the absolute isolation;
2. Remove the granulation tissue of the reabsorption area;
3. Neutralization of the medium with calcium hydroxide paste;
4. Remove the calcium hydroxide paste in the following session and execute endodontical treatment in the conventional way until the height of the reabsorption;
5. Apply BIO-C® REPAIR in the place of reabsorption, condensing it against the walls with specific instruments or a sterile ball of cotton slightly moistened;
6. Execute the coronary sealing with glass ionomer and restoration with the material preferred;
7. X-ray and execute preservation for at least two years;

D. TREATMENT OF EXTERNAL OR COMMUNICATING INTERNAL REABSORPTION SURGICALLY

1. Anesthetize, raise the gingival flap and execute osteotomy to locate the perforation;

2. Remove all the cement and dentin affected, removing the cells inducing reabsorption;
3. Execute hemostatic control and dry the cavity. It is advisable to condition the surface of the root affected, preferably with citric acid;
4. Insert BIO-C® REPAIR in the cavity which must be slightly moist; adapt the material with condensers and remove the excess if overflowing;
5. Reposition the gingival flap and sew it;
6. X-ray and execute preservation for at least two years.

Important: if there is external reabsorption above the bone crest it is not advisable to use BIO-C® REPAIR.

E. RETROFILLING IN PARAENDODONTIC SURGERY

1. Anesthetize, raise the gingival flap and execute osteotomy to expose the dental apex;
2. Cut 3 mm of the root apex which contains a large number of collateral canals (apical delta);
3. Execute retro preparation with appropriate instrument;
4. Execute hemostatic control and maintain the cavity slightly moist;
5. Insert BIO-C® REPAIR in the cavity, adapt with appropriate condensers and remove the excess overflowing material;
6. Further bleeding in the place from the periodontal ligament and bone tissue, return the gingival flap to its position and sew it;
7. X-ray and execute preservation for at least two years.

F. DIRECT AND INDIRECT PULP CAPPING

1. Anesthetize and install the absolute isolation;
 2. Remove the caries; if there is pulp exposure execute hemostasis;
 3. Apply BIO-C® REPAIR in the cavity or on the pulp, if there is exposure. Use appropriate instruments or a sterile ball of cotton slightly moistened for seating of the material;
 4. Recover BIO-C® REPAIR with a glass ionomer cement and wait 4 to 6 weeks to execute the final restoration. The glass ionomer can be used as the basis of the restoration;
- Warning: for Young permanent or deciduous teeth with a great exposure, execute the pulpotomy and then do steps 3 and 4 above.

G. APEXIFICATION

First Session

1. Anesthetize and install the absolute isolation;
2. Execute the preparation of the duct by the conventional technique;
3. Apply an intracanal dressing with calcium hydroxide paste for a week.

Second Session

1. Remove the intracanal dressing with files and irrigating solution;
2. Dry the duct with paper cones;
3. Make a mechanical barrier in the apical region to avoid overflowing of BIO-C® REPAIR. It is advisable to use a collagen sponge or similar products.
4. Apply BIO-C® REPAIR condensing it to form an apical cover of 3 to 4 mm;
5. X-ray to check the correct filling of the canal;
6. Place a sterile cotton ball slightly moistened at the entrance of the canal followed by provisional restoration.

Third Session

1. Remove the provisional restoration and fill the canal to the limit of the apical cover;
2. Restore with glass ionomer;
3. Execute the clinical and radiographic control of 3 to 6 months until the radiographic check of apical barrier of hard tissue.

H. APICIGENESIS AND PULPOTOMY

1. Anesthetize, install the absolute isolation and then make the access to the pulp chamber;
2. Remove the roof of the pulp chamber and execute excision of the coronary pulp with sharpened curettes or drills;
3. Cut the pulp 0.5 mm below the inlet of the root canal;
4. Wash the surgical wound abundantly with distilled water, physiological serum or specific detergent;
5. Execute hemostasis and maintain the cavity slightly moist;
6. Apply BIO-C® REPAIR on the pulp stump and adapt it with a sterile cotton ball slightly moistened;
7. X-ray to check the correct filling of the canal;
8. Restore the cavity provisionally with glass ionomer;
9. Execute the clinical and radiographic control of 3 to 6 months until the radiographic check of radicular formation;
10. After the end of the radicular formation choose between conventional endodontic treatment or only the restoration of the coronary endodontic cavity.

I. PULP REGENERATION

1. Anesthetize, install the absolute isolation and then make the access to the pulp chamber;
2. Empty the canal to the foramen, with the least wear possible of the target walls;
3. Irrigate abundantly with sodium hypochlorite at 2.5% and dry the canal with paper tips;
4. Place the Ca(OH)₂ paste in salt solution and execute the sealing with glass ionomer;

5. Maintain the intracanal medication for 4 weeks;
6. After 4 weeks, anesthetize the patient, install the absolute isolation, access the root canal and remove the paste $\text{Ca}(\text{OH})_2$ with irrigation of EDTA 17% (if possible, activate EDTA with Ultrasonic Tip in 3 cycles of 10 seconds' duration);
7. Take an endodontic file to the canal and further bleeding through injuring the tissue of the periapex to induce bleeding and later forming of coagulation;
8. Stable the bleeding 3 mm below the cemento-enamel junction and wait 15 minutes for the clot to form; if necessary place a collagen membrane at the start of the cervical third. This membrane will serve to contain the BIO-C® REPAIR in the cervical third;
9. Place a layer of 3 mm of BIO-C® REPAIR on the clot;
10. Place a sterile ball of cotton slightly moistened on the BIO-C® REPAIR, and wait 15 minutes for the start setting of the material;
11. Remove the ball of cotton, dry the cavity and apply a layer of glass ionomer. Then, restore the tooth with the material preferred professionally;
12. Execute the clinical and radiographic follow-up every 6 months until complete root formation.

WARNINGS

- Use absolute isolation in all the indications of user by canal;
- Close the syringe immediately after removing the material with appropriate pressure to avoid the hardening of the product;
- Before use check if the product is hardened at the tip of the syringe and, if so, discard the first portion;
- Never store the product in the refrigerator, as it is sensitive to humidity;
- The cavities where the product will be inserted do not need to be completely dry, but they must not be excessively moist;
- Do not use the product to fill the canals, as its flow is inappropriate for this purpose;
- Only use the product after controlling the infection, as in this phase the setting of the material does not occur due to the acid pH of the medium;
- Do not use the product when there is communication with the gingival groove, as the disintegration of the material will occur;
- Avoid causing excessive overflow of the product to the periapical region;
- Use sterile instruments to handle and apply the product;
- Product of single use;
- Clean and disinfect the syringe in ensuing uses to avoid cross contamination.

PRECAUTIONS

- Wear safety goggles, mask and gloves when handling the product. If there is contact with the eyes or skin, wash with water;
- Do not use the product after the validity term has expired;
- Do not use the product if the packaging is damaged;
- Do not use in patients with reported sensitivity to the components of the formula;
- Product not tested in pregnant and breast-feeding women;
- Keep it out of the reach of children.
- Before use check if the product is hardened at the tip of the syringe and, if so, discard the first portion.

STORAGE

- Maintain the product in a dry and ventilated place between 15°C and 30°C and with relative air humidity below 60%;
- Never store the product in a refrigerator, as it is sensitive to humidity;
- * Product of exclusive use of the Dental Surgeon.

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